

Joerns Custom Sling Order Form Oxford Deluxe Transport



STEP 3



Prior to requesting changes or if you're in any doubt, please consult this sling's standard specification on page 4.


STEP 2 Core Sling Dimensions: Select a standard size or enter custom dimensions.							
	Small	Medium	Large	Custom			
Dimension A	670mm	840mm	1100mm	mm			
Dimension B	290mm	340mm	400mm	mm			
Dimension C	860mm	1010mm	1060mm	mm			

Note: Use page 3 to illustrate bespoke sizing. Dimension A is the mid-point of the sling body. All dimensions are subject to manufacturing tolerances.

Material Options: Select the core material for your sling and include any bespoke requirements in the box provided.

Polyester Solid	Special Instructions:
Super-Soft	

Note: Some slings incorporate more than one fabric type - please consider this when requesting material changes. See page 4 for more information.



STEP 6 Head Supp	ort: Use page :	3 to help illustrate th	ne shape and size of	head support if required.	
Add Head Support	Add Adjustment Straps		Special Instruction	15:	
	Add Stiffer	ners			
]			
STEP 7 Straps, Loo	ops, Clips	& Handles:	Adjust strap lengths	number of loops, loop colours and positioning handles.	
Standard Strap Length	Custom Strap Length	Standard No. of Loops	No. of Extra Loops?	Special Instructions: (e.g. quantity and placement of positioning handles)	
Shoulder 510mm	mm	4			
Leg 540mm	mm	4			
Add Positioning Handle(s)		Clip Attachments	(Securi3®)		
Add Modesty Loop		When choosing clips, please specify			
		the length of strap required from the sling body to the clip (required for both shoulder and leg attachment points).			
Loop Colour Options: Black blue o	reen, grev orange			selection/sequence in the special instructions box above.	
See page 3 for dimensional illustration					
STEP 8 Body Belt:	Add or alter - L	what to anti-	upport and accounts	uring transfer. Use page 3 to help illustrate your needs.	
STEP 8 Body Beit.	Add or alter a boo				
Buckle Belt Length	Belt Wid	Ith Special I	nstructions: (e.g. pao	ding, elastic, etc)	
	1	nm			
STEP 9 Other Option	ONS: Use page	3 to help illustrate y	our bespoke require	ments.	
•			ecial Instructions:		
Add Positioning Pocket	Add Centre				
STEP 10 Additional	Requirem	ents			
STEP 11 Confirmation	on				
I am happy with and hereby gi	ve permission to p	roceed with the cus	stom sling order defir	ed on this form.	
Please call me to discuss/conf	firm this order and	or price prior to ma	nufacture. Pl	none No:	
L					
Customer PO Number:	Customer PO Number: Ordered By:				
Order Quantity:				Signed:	
Sales Order Number:				Date:	
Please Note: Joerns Healthcare reser	ve the right to refu	se a special order if	deemed unsafe. In si	uch event, we will contact you to discuss the matter further.	


Please return your completed form to:

orders@joerns.co.uk Please attach any supporting illustrations when submitting your form.

(T) +44 (0)344 811 1158 (W) www.joerns.co.uk

Illustrative Guidance

Please use this page to help illustrate or clarify any bespoke requirements



Standard Specification Breakdown

Please consider the following information when requesting changes to your standard sling.

	Sling Model			
	Oxford Deluxe Transport			
Material				
Core Material	Polyester Solid (Grey)			
Spine Overlay	Non-slip mesh overlay to inside centre of sling body			
Leg Overlay	N/A			
Padding				
Head Support	N/A			
Body	20mm Soft Foam + Hollow Fibre filling to under-arm pads			
Legs	10mm Soft Foam			
Head Support				
Included	N/A			
Integrated Stiffeners	N/A			
Adjustment Straps	N/A			
Loops, Clips & Handles				
Sling Attachment Type	Loops			
Positioning Handles	N/A			
Modesty Loop	N/A			
Body Belt				
Included	\checkmark			
Attachment Type	Buckle			
Other Options				
Centre Markers	N/A			
Positioning Pocket	N/A			

You do not need to submit this page, it is for your reference only.